

Brava[®] Skin Barrier


An invisible friend for the skin



Shields against output and adhesives

Leakage and adhesives can cause a variety of skin problems. Brava Skin Barrier leaves a thin breathable film on the skin that protects against potential damage without affecting the application of the next appliance.

The Brava Skin Barrier is a silicone-based product so it does not sting when applied. It comes as a spray and in travel-friendly wipes.



Now with
even more
liquid



Brava[®] Skin Barrier

Prevents residue build-up


Brava Skin Barrier reduces the risk of residue when removing adhesive. It either comes off when removing the adhesive or it can be removed by gentle cleaning. This means residue does not build up on the skin over time.

Dries in seconds

With Brava Skin Barrier, your patients will be ready to apply a new adhesive almost immediately.

Product name	Item number	Unit in package	Volume
Brava Skin Barrier spray	12020	1	50ml
Brava Skin Barrier wipe	12021	30	-



 Before use, always consult the "instructions for use" document delivered with the products.

Brava supporting products

For customising appliance fit and for issues such as leakage and skin irritation.

What's right for your patients?



Brava[®]

Ostomy Care / Contenance Care / Wound & Skin Care / Urology Care

Coloplast A/S, Høtveddam 1, 3050 Humlebaek, Denmark
www.coloplast.com The Coloplast logo is a registered trademark of Coloplast A/S. © [YYYY-MM.] All rights reserved Coloplast A/S. PM-03910

